

Bond strength of universal adhesive applied to dry and wet dentin: one-year in vitro evaluation

Maria Luísa de Alencar e Silva Leite¹, Caroline de Farias Charamba², Renally Bezerra Wanderley e Lima³, Sônia Saeger Meireles², Rosângela Marques Duarte², Ana Karina Maciel de Andrade^{2,*}

¹ Department of Dental Materials and Prosthodontics, Universidade Estadual de São Paulo (Unesp), School of Dentistry, Araraquara, SP, Brazil.

² Department of Restorative Dentistry, Universidade Federal da Paraíba, João Pessoa, PB, Brazil.

³ Department of Restorative Dentistry, Piracicaba Dental School, Universidade de Campinas (UNICAMP), Piracicaba, São Paulo, Brazil

Corresponding author:

Ana Karina Maciel de Andrade
 <https://orcid.org/0000-0003-4520-5176>
Department of Restorative Dentistry, Federal University of Paraíba, João Pessoa, PB, Brazil. Cidade Universitária, s/n, Castelo Branco, João Pessoa, PB, Brazil. Postal code: 58051-900. Phone: +55 83 999411191. E-mail: kamandrade@hotmail.com

Received: May 05, 2019

Accepted: October 10, 2019

Aim: This study evaluated the influence of dentin wettability on the immediate and extended microtensile bond strength (mTBS) of a universal adhesive system used in the etch-and-rinse strategy. **Methods:** Twenty human third molars were selected and divided into four groups according to the adhesive system and dentin wettability. The mTBS values of each group were registered 24 h and one year after adhesive system application and resin composite block build-up (n=30). Data were analyzed by the t-test ($p < 0.05$). **Results:** When both adhesive systems were compared, there was no statistically significant difference when they were applied following wet bonding ($p > 0.05$). However, the dry bonding reduced μ TBS values of the Adper Single Bond 2 adhesive ($p < 0.05$). Regarding storage time, both groups presented similar μ TBS values at 24 h and one year ($p > 0.05$). **Conclusions:** Therefore, the Scotchbond Universal Adhesive can be applied to dry or wet dentin without compromising the etch-and-rinse bonding quality and the durability of the restorations.

Keywords: Dental cements. Dentin. Tensile strength.



Introduction

Due to the high patient demand for esthetic restorations and minimally invasive treatments, adhesive systems have become a fundamental material for achieving the stable and long-term bonding effectiveness of aesthetic restorations to mineralized tooth tissues¹. Adhesion of restorative materials to hard tissues (enamel and dentin) is challenging due to the differences in these tissues². Dental enamel presents a structure more homogenous than that of dentin, resulting in reliable and long-term durable adhesive bonds between the enamel and the restorative material. In contrast, dentin substrate is characterized by a wide variety of inorganic and organic components; consequently, adhesive bonding to dentin is more sensitive².

In the past, the adhesive systems available on the market have been classified into two categories: etch-and-rinse (ER) and self-etch strategies (SE)^{3,4}. In the ER strategy, a phosphoric acid gel is first applied to the dentin substrate, followed by application of the primer and the bond resin separately or in a single solution^{3,4}. A major disadvantage of this strategy is susceptibility to variations in the degree of dentin moisture, which is subjective and depends on operator skills⁵. Excessive residual moisture may hinder the impregnation of monomers into demineralized substrate by dilution of these components⁶. Conversely, the overdrying of the dentin surface may promote collagen fibril collapse and, consequently, the incomplete impregnation of resin monomers into the collagen fibers, decreasing the bond strength⁷. In addition, the dentin etching with phosphoric acid can lead to collagen matrix degradation in the dentin as a result of the activation of endogenous dentin collagenolytic enzymes by acidity^{8,9}. This process can result in impaired bond integrity of the adhesive interface^{8,9}.

To overcome the limitations of the ER strategy, SE adhesives were developed, characterized by acidic functional monomers that simultaneously etch and prime the tooth substrate for bonding¹⁰. Previous studies have demonstrated that the laboratory and clinical performance of traditional SE adhesives has not been satisfactory, mainly in the enamel^{11,12}. Selective enamel etching has been suggested to improve the adhesive bond to enamel^{13,14}. However, the application of phosphoric acid to enamel may inadvertently etch dentin, decreasing the bond strength of SE adhesives to this tissue¹⁵.

The latest trends in adhesive systems are multi-mode or universal one-bottle adhesives that may be used as either ER or SE adhesives or as SE adhesives on dentin and ER adhesives on enamel (selective enamel-etching)¹⁶⁻¹⁸. Most of the multi-mode or universal adhesives are characterized by 10-methacryloyloxydecyl dihydrogen phosphate (10-MDP), which bonds ionically to dentin, forming hydrolytically stable calcium salts on hydroxyapatite in the form of "nano-layering"¹⁶. The manufacturers claim that multi-mode adhesives may also be applied to dentin under different bonding strategies (dry or wet bonding)^{15,19}. However, there are few studies of the bonding performance and reliability of those adhesives over the short or long term when applied to dry or wet dentin^{15,19,20}. Thus, more information is necessary to predict the long-term bonding durability of universal adhesives.

Therefore, the aim of this *in vitro* study was to evaluate whether the bond strength of a universal adhesive to dentin is affected by application mode (dry- and wet-bonding),

when the etch-and-rinse strategy is used, after 24 h and one year of water storage. Three null hypotheses were set: 1) Bond strength to dentin would not be affected by application mode (dry- and wet-bonding) when the etch-and-rinse protocol is used; 2) bond strength to dentin would not be affected by storage period; and 3) there would be no difference in bond strength to dentin between etch-and-rinse and universal adhesives.

Materials and methods

Tooth selection and preparation

Twenty healthy human third molars were used in this study, after approval from the Research Ethics Committee of the University of Paraiba (protocol n. 17665613.2.0000.5188). The teeth were cleaned, stored in 0.2% thymol solution, and used within one month after extraction. All tooth roots were embedded in self-curing acrylic resin. Then, the occlusal enamel was removed by means of a diamond disc (Extec, Enfield, CT, USA) under water-cooling. The exposed occlusal dentin surfaces were wet-abraded with silicon carbide paper (600 grit) under water-cooling for 60 s by means of a polishing machine (Politriz ERIOS – 27000, São Paulo, SP, Brasil) to standardize the smear layer²¹.

Experimental design

The teeth were randomly assigned among four groups according to the different bonding strategies of the selected adhesive systems (n = 5). The two-step etch-and-rinse adhesive, Adper Single Bond 2 (AS) (3M ESPE, St. Paul, MN, USA), and the universal adhesive, Scotchbond Universal Adhesive (SU) (3M ESPE, St. Paul, MN, USA), were applied to dentin surfaces following a dry- or a wet-bonding etch-and-rinse adhesive protocol. Composition, batch number of each material, and adhesive strategies are shown in Table 1.

Table 1. Brand, batch number, composition, and adhesive strategies of materials used.

Adhesive system	Manufacturer/ batch number	Type	Composition	Etch-and-rinse strategy
Scotchbond Universal	3M/ESPE, St. Paul, MN, USA (526247)	Universal adhesive system	10-MDP, phosphate monomer, dimethacrylate resins, HEMA, methacrylate-modified polyalkenoic acid copolymer, filler, ethanol, water, initiators, silane	I. Apply etchant for 15 s. II. Rinse for 10 s. III. Air-dry to remove excess water. IV. Keep dentin moist (wet-bonding approach) or keep dentin dry. Do not overdry (dry-bonding approach). V. Apply 2 consecutive coats of adhesive. VI. Gently air-dry for 5 s. VII. Light-polymerize for 10 s.
Adper Single Bond 2	3M/ESPE, St. Paul, MN, USA (N49344)	Etch-and-rinse adhesive system	1. Etchant: 35% phosphoric acid (Scotchbond Etchant). 2. Adhesive: bis-GMA, HEMA, dimethacrylates, ethanol, water, photoinitiator, methacrylate functional copolymer of polyacrylic and poly(itaconic) acids, 10% by weight of 5- nm-diameter spherical silica particles	I. Apply etchant for 15 s. II. Rinse for 10 s. III. Air-dry to remove excess water. IV. Keep dentin moist (wet-bonding approach) or keep dentin dry. Do not overdry (dry-bonding approach). V. Apply 2 consecutive coats of adhesive. VI. Gently air-dry for 5 s. VII. Light-polymerize for 10 s.

Restorative procedure and specimen preparation

After the bonding process, three resin composite increments (Z100-3M ESPE, St. Paul, MN, USA) of 1.5 mm were placed on the dentin surface, and each increment was light-cured for 40 s by means of a LED light-curing unit set at 400 mW/cm² (GNATUS, São Paulo, SP, Brazil). The restored teeth were then stored in distilled water at 37°C ($\pm 1^\circ\text{C}$) for 24 h²².

After this storage period, the specimens were sectioned longitudinally in the mesio-distal and buccal-lingual directions across the bonded interface, by means of a slow-speed diamond disc (Labcut 1010, Extec, Enfield, CT, USA). Then, the specimens were sectioned transversely in the cervical region to obtain bars measuring 1 mm² x 10 mm²¹. Half of the bars obtained from each tooth were used immediately (24 h) for the micro-tensile bond strength (μTBS) test, while the other half were stored in distilled water at 37°C ($\pm 1^\circ\text{C}$) for one year²² and then subjected to μTBS testing. For μTBS testing, the bars were fixed to a testing jig with cyanoacrylate glue (Super Bonder Gel – Loctite Brasil Ltda) and subjected to tensile load at a crosshead speed of 0.5 mm/min until failure (microtensor OM-100 machine, Odeme, Luzerna, SC, Brasil)²³. The μTBS values (MPa) were calculated by dividing the load at failure by the cross-sectional bonding area. The fractured surfaces of all specimens were observed by means of an optical microscope (XJM-400, KOZO, Nanjing, China) at a magnification of 100x, and fracture patterns were classified as (1) cohesive failure in adhesive, (2) cohesive failure in dentin, (3) cohesive failure in the hybrid layer, or (4) mixed failure (cohesive failure in adhesive and cohesive failure in the hybrid layer). The data from fracture patterns were analyzed by descriptive statistics.

The experimental unit in the current study was the bar. An average of 10 to 15 bars was obtained from each tooth, with the experimental group ($n = 30$) having the smallest number of test specimens, the experimental groups are in the figure 1. Thus, the μTBS value of each bar was used for statistical analysis. The data from μTBS were analyzed by the t -test for independent samples ($\alpha = 0.05$). In addition, the reliability of the bond




Figure 1. Sample diagram.

strength for each group was analyzed by Weibull analysis. The Weibull moduli (shape parameter) (slope of the line relating applied stress and the probability of specimen failure, m) were calculated, applying maximum likelihood estimation. The 95% upper and lower confidence intervals were calculated using the likelihood ratio (MINITAB 17.0, State College, Pennsylvania, USA). Differences between the paired values for m were considered significant when the 95% confidence intervals did not overlap.

Results

The bond strength means (MPa) and standard deviations for all experimental groups are presented in Table 2.

Table 2. Microtensile bond strength (MPa) values of adhesive systems to dentin among the test groups, comparing materials and different etch-and-rinse strategies.

Storage time	Etch-and-rinse strategy	Adhesive system	
		Scotchbond Universal (SU)	Adper Single Bond 2 (AS)
24 hours	wet-bonding	49.08 (15.23) ^{Aa}	47.38 (23.23) ^{Aa}
	dry-bonding	52.28 (19.43) ^{Aa}	26.91 (10.25) ^{Bb}
One year	wet-bonding	50.22 (16.36) ^{Aa}	45.13 (18.61) ^{Aa}
	dry-bonding	53.87 (23.51) ^{Aa}	21.63 (13.91) ^{Bb}

Means followed by different uppercase letters in the same row indicate statistically significant differences between adhesive systems within the same etch-and-rinse strategy and evaluation time. Means followed by different lowercase letters on the same column indicate statistically significant differences between etch-and-rinse strategy within the same adhesive system and evaluation time ($p < 0.05$).

The SU adhesive showed no difference in μ TBS values between dry- and wet-bonding etch-and-rinse strategies in 24 hours and after one year ($p=0.48$). Conversely, the AS adhesive applied to a dry dentin surface presented significantly lower μ TBS values than that measured when the adhesive was applied to the dentin following a wet-bonding technique ($p < 0.001$). When both adhesive systems were compared, there was no statistically significant difference between them when they were applied following a wet-bonding strategy in 24 hours ($p=0.74$) and after one year ($p=0.26$). However, when the adhesive systems were applied to dry dentin, the SU showed significantly higher μ TBS values ($p < 0.001$) for both storage times.

Regarding storage time, both adhesive systems presented similar μ TBS values at 24 hours and one year ($p > 0.05$), regardless of the etch-and-rinse strategy (dry or wet dentin).

The Weibull analysis indicated no change in the Weibull parameter for SU adhesive for all study conditions. On the other hand, AS adhesive showed the lowest m values applied to the dentin following a wet-bonding (24 hours) and dry-bonding (one year) technique. Comparing the adhesive systems, SU demonstrated higher m values than

AS for wet-bonding strategy in 24 hours and dry-bonding strategy in one year. This finding suggests a more predictable, consistent performance of SU product than seen when using the AS adhesive (Table 3, Figure 2, and 3).

Table 3. Weibull moduli (m) values, among the experimental groups comparing different adhesive system, storage time, and etch-and-rinse strategy.

Storage time	Etch-and-rinse strategy	Adhesive system	
		Scotchbond Universal (SU)	Adper Single Bond 2 (AS)
24 hours	wet-bonding	3.6 (2.73-4.74) ^{Aa}	2.23 (1.67-2.98) ^{Bab}
	dry-bonding	3.12 (2.32-4.2) ^{Aa}	2.80 (2.16-3.63) ^{Aa}
One year	wet-bonding	3.37 (2.57-4.42) ^{Aa}	2.68 (2.02-3.56) ^{Aa}
	dry-bonding	2.57 (1.9-3.46) ^{Aa}	1.68 (1.28-2.19) ^{Bb}

Means followed by different uppercase letters in the same row indicate statistically significant differences between adhesive systems within the same etch-and-rinse strategy and evaluation time. Means followed by different lowercase letters on the same column indicate statistically significant differences between etch-and-rinse strategy within the same adhesive system and evaluation time ($p < 0.05$).




Figure 2. Weibull distribution plots of microtensile bond strength data for different etch-and-rinse strategy and evaluation time within the same adhesive system. ASW- Adper Single Bond 2 wet-bonding; ASD- Adper Single Bond 2 dry-bonding.

The percentages of specimens according to fracture mode for all experimental groups at 24 h and one year are summarized in Tables 4 and 5, respectively. The predominant failure patterns observed were mixed failure, cohesive failure in adhesive, and cohesive failure in the hybrid layer (IV).




Figure 3. Weibull distribution plots of microtensile bond strength data for different etch-and-rinse strategy and evaluation time for Scotchbond Universal adhesive system. SUW- Scotchbond Universal wet-bonding; SUD- Scotchbond Universal dry-bonding.

Table 4. Percentages (%) of specimens according to the fracture mode of all test groups in 24 hours' storage time.

Adhesive systems	Etch-and-rinse strategy	Fracture patterns (%)			
		I	II	III	IV
Scotchbond Universal	wet-bonding	23.3	0	3.3	73.4
	dry-bonding	13.3	0	13.3	73.4
Adper Single Bond 2	wet-bonding	16.5	0	16.25	66.7
	dry-bonding	13.3	0	20	66.7

Table 5. Percentages (%) of specimens according to the fracture mode of all test groups in one-year storage time.

Adhesive systems	Etch-and-rinse strategy	Fracture patterns (%)			
		I	II	III	IV
Scotchbond Universal	wet-bonding	0	0	0	100
	dry-bonding	20	0	0	80
Adper Single Bond 2	wet-bonding	13.3	0	0	86.7
	dry-bonding	50	0	0	50

Discussion

Based on the results of this in vitro study, the universal adhesive etch-and-rinse strategy can be applied to dry and wet dentin surfaces, since dry- and wet-bonding strategies demonstrated suitable bond strength to dentin in short (24 h) and long periods of aging (one year). In contrast, the dry-bonding mode had a negative impact on dentin bond quality of the two-step adhesive, and, consequently, this material presented the worst bond strength values when applied to dry dentin. Thus, the null hypotheses of the present study were rejected.

The universal adhesives were introduced to simplify and optimize the adhesive procedures. These materials can be used in different types of adherent substrate and etching modes (etch-and-rinse or self-etch). A recent systematic review and meta-analysis showed revealed that the mild universal adhesives seem stable materials, in both etch-and-rinse or self-etch strategies. In this study, the universal adhesive was used in etch-and-rinse mode. As is known, the SU adhesive tested in this study contains 10-MDP monomer, which provides chemical bonding to hydroxyapatite, forming hydrolytically stable calcium salts. Additionally, this adhesive contains a polyalkenoic acid copolymer that also interacts ionically with the hydroxyapatite through the carboxylic groups^{24,25}. To trigger the ionization of the phosphate monomer and polyalkenoic acid copolymer, a certain amount of water (10–15% by wt) (3M ESPE) was added to the universal adhesive²⁶. Thus, the unique composition of the universal adhesive (SU) can explain the results observed in this study, which revealed no statistically significant difference for μ TBS values between the dry- or wet-bonding etch-and-rinse strategy after 24 hours and one year. The chemical components and the water contained in SU adhesives may be able to rehydrate collagen fibrils, allowing for the re-expansion of the interfibrillar spaces for the infiltration of resin monomers and the formation of resin tags branching out profusely into dentin tubules^{27,28}, permitting a satisfactory dentin seal and performance of those adhesives, even with dry-demineralized dentin^{15,19,29}.

The results of this study are in agreement with those of other reports^{15,19,25,28-31}, which concluded that the bonding durability of the universal adhesive was acceptable and did not seem to vary depending on the etching and application mode (dry- or wet-bonding). This outcome was confirmed by Weibull analysis that revealed how reliable a given treatment is³². The high m values observed for SU adhesive (table 3, figure 1), mainly with dry-demineralized dentin, suggest that the universal adhesive investigated can be considered reliable over time.

Conversely, μ TBS values for the conventional two-step etch-and-rinse adhesive (AS) decreased when this adhesive was applied to a dry-dentin surface in the etch-and-rinse protocol after 24 h and one year. This can be attributed to an insufficient amount of water within the dentin structure, which can lead to a collapsed collagen network⁷. Consequently, reduced resin monomer penetration into the entire depth of the decalcified dentin can occur when the AS adhesive is used as the etch-and-rinse adhesive to dry dentin, corroborating with other studies⁷. The low m values observed for the AS adhesive with dry dentin (table 3, figure 2) ratify that the application of this material to wet dentin is recommended to obtain mechanical stability of the bond to dentin.

Regarding failure patterns, mixed failure (cohesive failure in adhesive and cohesive failure in the hybrid layer) was the predominant mode for all experimental groups. This finding is in agreement with other studies^{5,33-36} and suggests that the bonding agents interacted with dental substrates during monomer infiltration by forming a hybrid layer. However, this layer fractured due to concentrated tension at the adhesive interface.

One of the main degradation mechanisms of dentin bonding is the hydrolysis of the collagen fibrils and the polymerized resin matrix in the adhesive layer³⁷. However, the patterns of degradation of resin-dentin bonding depend on the type of adhesive system³⁷. In this study, μ TBS values for both adhesives tested did not decrease significantly after one year of storage in distilled water. Therefore, the storage medium used in this research did not influence the long-term stability of the bond-adhesive interface.

Therefore, according to the present study, the SU adhesive may be an effective alternative approach for dentin restorations. Nevertheless, other factors related to the oral environment may have a more complex influence on the bonding performance of adhesives. Thus, future investigations should be conducted for the clinical evaluation of the SU adhesive to confirm the results presented in this research.

In conclusion, the universal adhesive tested in this study can be applied to dry or wet demineralized dentin without compromising the etch-and-rinse bonding quality and the durability of the restorations. On the other hand, the dry-bonding etch-and-rinse adhesive protocol influenced negatively the dentin bond quality of the conventional two-step adhesive.

Acknowledgments

The authors would like to acknowledge the financial support from the Conselho Nacional de Desenvolvimento Científico e Tecnológico (CNPq, Brazil).

References

1. Takamizawa T, Barkmeier WW, Tsujimoto A, Berry TP, Watanabe H, Erickson RL, et al. Influence of different etching modes on bond strength and fatigue strength to dentin using universal adhesive systems. *Dent Mater*. 2016 Feb;32(2):e9-21. doi: 10.1016/j.dental.2015.11.005.
2. Haller B. Recent developments in dentin bonding. *Am J Dent*. 2000 Feb;13(1):44-50.
3. De Munck J, Van Den Steen PE, Mine A, Van Landuyt KL, Poitevin A, Opendakker G, et al. Inhibition of enzymatic degradation of adhesive-dentin interfaces. *J Dent Res*. 2009 Dec;88(12):1101-6. doi: 10.1177/0022034509346952.
4. Reis A, Pellizzaro A, Dal-Bianco K, Gones OM, Gones OM, Patzlauff R, Loguercio AD. Impact of adhesive application to wet and dry dentin on long-term resin-dentin bond strengths. *Oper Dent*. 2007 Jul-Aug;32(4):380-7.
5. Tay FR, Pashley DH. Aggressiveness of contemporary self-etching systems. I: Depth of penetration beyond dentin smear layers. *Dent Mater*. 2001 Jul;17(4):296-308.
6. Van Meerbeek B, De Munck J, Mattar D, Van Landuyt K, Lambrechts P. Microtensile bond strengths of an etch&rinse and self-etch adhesive to enamel and dentin as a function of surface treatment. *Oper Dent*. 2003 Sep-Oct;28(5):647-60.

7. Da Silva MA, Rangel PM, Barcellos DC, Pagani C, Rocha Gomes Torres C. Bond strength of adhesive systems with different solvents to dry and wet dentin. *J Contemp Dent Pract.* 2013 Jan 1;14(1):9-13.
8. Mazzoni A, Scaffa P, Carrilho M, Tjäderhane L, Di Lenarda R, Polimeni A. Effects of etch-and-rinse and self-etch adhesives on dentin MMP-2 and MMP-9. *J Dent Res.* 2013 Jan;92(1):82-6. doi: 10.1177/0022034512467034.
9. Tjäderhane L, Nascimento FD, Breschi L, Mazzoni A, Tersariol IL, Geraldini S, et al. Strategies to prevent hydrolytic degradation of the hybrid layer—a review. *Dent Mater.* 2013 Oct;29(10):999-1011. doi: 10.1016/j.dental.2013.07.016.
10. Salz U, Bock T. Testing adhesion of direct restoratives to dental hard tissue—a review. *J Adhes Dent.* 2010 Oct;12(5):343-71. doi: 10.3290/j.jad.a19741.
11. Peumans M, Kanumilli P, De Munck J, Van Landuyt K, Lambrechts P, Van Meerbeek B. Clinical effectiveness of contemporary adhesives: a systematic review of current clinical trials. *Dent Mater.* 2005 Sep;21(9):864-81.
12. Perdigão J, Carmo AR, Anauate-Netto C, Amore R, Lewgoy HR, Cordeiro HJ, et al. Clinical performance of a self-etching adhesive at 18 months. *Am J Dent.* 2005 Apr;18(2):135-40.
13. Frankenberger R, Lohbauer U, Roggendorf MJ, Naumann M, Taschner M. Selective enamel etching reconsidered: better than etch-and-rinse and self-etch? *J Adhes Dent.* 2008 Oct;10(5):339-44.
14. Jacker-Guhr S, Sander J, Luehrs AK. How “Universal” is Adhesion? Shear Bond Strength of Multi-mode Adhesives to Enamel and Dentin. *J Adhes Dent.* 2019;21(1):87-95. doi: 10.3290/j.jad.a41974.
15. Mena-Serrano A, Kose C, De Paula EA, Tay LY, Reis A, Loguercio AD, et al. A new universal simplified adhesive: 6-month clinical evaluation. *J Esthet Restor Dent.* 2013 Feb;25(1):55-69. doi: 10.1111/jerd.12005.
16. Yoshida Y, Yoshihara K, Nagaoka N, Hayakawa S, Torii Y, Ogawa T, et al. Self-assembled nano-layering at the adhesive interface. *J Dent Res.* 2012 Apr;91(4):376-81. doi: 10.1177/0022034512437375.
17. Lenzi TL, Soares FZ, Raggio DP, Pereira GK, Rocha RO. Dry-bonding etch-and-rinse strategy improves bond longevity of a universal adhesive to sound and artificially-induced caries-affected primary dentin. *J Adhes Dent.* 2016;18(6):475-82. doi: 10.3290/j.jad.a36670.
18. Cuevas-Suárez CE, da Rosa WLO, Lund RG, da Silva AF, Piva E. Bonding performance of universal adhesives: an updated systematic review and meta-analysis. *J Adhes Dent.* 2019;21(1):7-26. doi: 10.3290/j.jad.a41975.
19. Hanabusa M, Mine A, Kuboki T, Momoi Y, Van Ende A, Van Meerbeek B, et al. Bonding effectiveness of a new ‘multi-mode’ adhesive to enamel and dentine. *J Dent.* 2012 Jun;40(6):475-84. doi: 10.1016/j.jdent.2012.02.012.
20. Perdigão J, Kose C, Mena-Serrano AP, De Paula EA, Tay LY, Reis A, et al. A new universal simplified adhesive: 18-month clinical evolution. *Oper Dent.* 2014 Mar-Apr;39(2):113-27. doi: 10.2341/13-045-C.
21. Orellana N, Ramírez R, Roig M, Giner L, Mercade M, Durán F, et al. Comparative study of the microtensile bond strength of three different total etch adhesives with different solvents to wet and dry dentin. (in vitro test). *Acta Odontol Latinoam.* 2009;22(1):47-56.
22. Shinohara MS, De Goes MF, Schneider LF, Ferracane JL, Pereira PN, Di Hipólito V, et al. Fluoride-containing adhesive: Durability on dentin bonding. *Dent Mater.* 2009 Nov;25(11):1383-91. doi: 10.1016/j.dental.2009.06.011.
23. Di Hipólito V, Reis AF, Mitra SB, de Goes MF. Interaction morphology and bond strength of nanofilled simplified-step adhesives to acid etched dentin. *Eur J Dent.* 2012 Oct;6(4):349-60.

24. Yoshida Y, Van meerbeek B, Nakayama Y, Snauwaert J, Hellemans L, Lambrechts P, et al. Evidence of chemical bonding at biomaterial-hard tissue interface. *J Dent Res.* 2000 Feb;79(2):709-14.
25. Vermelho PM, Reis AF, Ambrosano GMB, Giannini M. Adhesion of multimode adhesives to enamel and dentin after one year of water storage. *Clin Oral Investig.* 2017 Jun;21(5):1707-5. doi: 10.1007/s00784-016-1966-1.
26. Luque-Martinez IV, Perdigão J, Muñoz MA, Sezinando A, Reis A, Loguercio AD. Effects of solvent evaporation time on immediate adhesive properties of universal adhesives to dentin. *Dent Mater.* 2014 Oct;30(10):1126-35. doi: 10.1016/j.dental.2014.07.002.
27. Dal-Bianco K, Pellizzaro A, Patzlaft R, de Oliveira Bauer JR, Loguercio AD, Reis A. Effects of moisture degree and rubbing action on the immediate resin-dentin bond strength. *Dent Mater.* 2006 Dec;22(12):1150-6.
28. Lenzi TL, Soares FZM, De Oliveira RR. Does Bonding Approach Influence the Bond Strength of Universal Adhesive to Dentin of Primary Teeth? *J Clin Pediatr Dent.* 2017;41(3):214-8. doi: 10.17796/1053-4628-41.3.214.
29. Loguercio AD, De Paula EA, Hass V, Luque-Martinez I, Reis A, Perdigão J. A new universal simplified adhesive: 36-Month randomized double-blind clinical trial. *J Dent.* 2015 Sep;43(9):1083-92. doi: 10.1016/j.jdent.2015.07.005.
30. Tsujimoto A, Shimatani Y, Nojiri K, Barkmeier WW, Markham MD, Takamizawa T, et al. Influence of surface wetness on bonding effectiveness of universal adhesives in etch-and-rinse mode. *Eur J Oral Sci.* 2019 Apr;127(2):162-9. doi: 10.1111/eos.12596.
31. Sugimura R, Tsujimoto A, Hosoya Y, Fischer NG, Barkmeier WW, Takamizawa T, et al. Surface moisture influence on etch-and-rinse universal adhesive bonding. *Am J Dent.* 2019 Feb;32(1):33-8.
32. Weibull W. A statistical distribution function of wide applicability. *J Appl Mech.* 1951 Set;18:293-7.
33. Hass V, Dobrovolski M, Zander-Grande C, Martins GC, Gordillo LA, Rodrigues Accorinte Mde L, et al. Correlation between degree of conversion, resin–dentin bond strength and nanoleakage of simplified etch-and-rinse adhesives. *Dent Mater.* 2013 Sep;29(9):921-8. doi: 10.1016/j.dental.2013.05.001.
34. Marchesi G, Frassetto A, Mazzoni A, Apolonio F, Diolosà M, Cadenaro M, et al. Adhesive performance of a multi-mode adhesive system: 1-year in vitro study. *J Dent.* 2014 May;42(5):603-12. doi: 10.1016/j.jdent.2013.12.008.
35. Munoz MA, Luque I, Hass V, Reis A, Loguercio AD, Bombarda NH. Immediate bonding properties of universal adhesives to dentine. *J Dent.* 2013 May;41(5):404-11. doi: 10.1016/j.jdent.2013.03.001.
36. Takai T, Hosaka K, Kambara K, Thitthaweerat S, Matsui N, Takahashi M, et al. Effect of air-drying dentin surfaces on dentin bond strength of a solvent-free one-step adhesive. *Dent Mater J.* 2012;31(4):558-63.
37. Breschi L, Mazzoni A, Ruggeri A, Cadenaro M, Di Lenarda R, de Stefano Dorigo E. Dental adhesion review: aging and stability of the bonded interface. *Dent Mater.* 2008 jan;24(1):90-101.